

PROCEDURE FOR USE OF AED

- | |
|---|
| ➤ If someone collapses – tap the victim and shout “Are you OK?” |
| ➤ If no response & the victim is not breathing or not breathing normally |
| ➤ Call 9-1-1 and get the AED |
| ➤ If someone else is present, send them to call 9-1-1 and get the AED |
| ➤ Begin Hands Only CPR – Push hard and fast in the center of the chest |
| ➤ Use the AED as soon as it arrives <ul style="list-style-type: none">○ Open the lid and turn on the AED○ Follow the voice prompts<ul style="list-style-type: none">➤ Expose the patient’s chest –if wet, wipe dry➤ Apply the electrode pads to the patient’s exposed chest, as shown in the pictures➤ Listen to the voice prompts –Clear (do not touch) the patient when instructed to do so➤ Press the SHOCK button, if instructed to do so |
| ➤ Restart compressions on chest |
| ➤ Continue to follow the voice prompts |
| ➤ If the patient starts moving or breathing normally, stop compressions and place the patient in the recovery position (on left side) |

KEEP AED ATTACHED UNTIL EMS ARRIVES